

Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

Immunoassay Batch Review Checklist

Batch Name: <u>EDA-20180626U-AAJ</u>		Analyst Review	Technical Review
Worklist	Batch name written and item(s) tested for each case listed	✓	✓
	Pipette(s) used listed	✓	✓
	Lot numbers for calibrators, controls, and reagents are listed	✓	✓
	Inventory number(s) listed	✓	✓
Sequence	Verify batch file name is consistent	✓	✓
	Reviewer verified, initialed, and dated sequence	✓	✓
	Verify appropriate calibrator and controls were analyzed, and appropriate controls were run at least after every 10 case samples	✓	✓
	Verify all cases were processed and have data	N/A	N/A
Magellan Data	Verify average absorbance of blank replicates is above 1.000	✓ <i>avg 1.2108</i>	✓
	Verify calibrators and controls are acceptable	✓	✓
	Verify case samples are bracketed by acceptable positive controls	N/A	N/A
	Verify correct method was run	✓ <i>avg 1.2108</i>	✓
Immunoassay Case Summary Report	Verify absorbance readings on Immunoassay Case Summary Report against Magellan data	N/A	N/A
	Verify results on Immunoassay Case Summary Report against Magellan data	↓	↓
	Verify analyst, matrix, instrument, and batch name(s)	↓	↓
	Verify forensic case number and item tested	↓	↓
Immunoassay Batch QC Data Worksheet	Verify batch and analyst names	✓	✓
	Verify preparation of calibrators, controls, and reagents	✓	✓
	Verify pipette(s) against those listed on the worklist	✓	✓
	Verify lot numbers and expiration dates of calibrators, controls, and reagents	✓	✓
	Verify kit lot numbers and expirations dates against Drug Standard Database	✓	✓
Verify all quality controls were added to QC Logs		✓	✓
All comments and/or strikethroughs, if any, initialed		✓	✓
All pages contain batch name and are initialed		✓	✓

AAJ

Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

Notes from Technical Review

CIA 20180626U-AAJ


- Check marks left off checklist "Magellan Data" section.
- Incorrect batch name on QC data sheet.
- Values transcribed incorrectly in QC log: carisoprodol (all), methadone (neg abs.), and oxycodone (cal to binding).

- emc 06/27/18


Analyst Review

Technical Review

Signature/Date:

 06/26/18

Signature/Date:

 06/27/18



Houston Forensic Science Center
Comparative and Analytical Division - Toxicology

Immunoassay Batch QC Data

Batch: EIA-2018 062601
~~EIA 06268~~ - AAS Analyst: Ashley Ann Johnson
ay 6/27/18

Matrix: Blood Urine

Instrument: Tecan Tecan-2

Pipette(s): 0395 235D

Phosphate Buffer Saline Lot Number: GA0918-PBS7 Expiration Date: 10/9/18

Controls:	Lot Number	Expiration	Oxycodone:	Lot Number	Expiration
Blank	<u>101017U</u>	<u>10/10/18</u>	Negative	<u>050718U-Q-50</u>	<u>10/9/18</u>
Negative	<u>050718U-MKH</u>	<u>10/27/18</u>	Calibrator	<u>050418U-C-100</u>	<u>10/9/18</u>
Calibrator	<u>050418U-MXC</u>	<u>10/27/18</u>	Positive	<u>050718U-Q-200</u>	<u>10/9/18</u>
Positive	<u>050718U-MXP</u>	<u>10/27/18</u>			

	Kit Lot Number	Expiration		Kit Lot Number	Expiration
Amphetamine	<u>EK 16566</u>	<u>2/19</u>	Methadone	<u>EK 16418</u>	<u>12/18</u>
Barbiturates	<u>EK 16588</u>	<u>2/19</u>	Methamphetamine	<u>EK 16296</u>	<u>11/18</u>
Benzodiazepines	<u>EK 16701</u>	<u>3/19</u>	Opiates	<u>EK 16544</u>	<u>2/19</u>
Benzoylcegonine	<u>EK 16643</u>	<u>3/19</u>	Phencyclidine	<u>EK 16593</u>	<u>2/19</u>
Cannabinoids	<u>EK 16653</u>	<u>3/19</u>	Zolpidem	<u>N/A</u>	<u>N/A</u>
Carisoprodol	<u>EK 16763</u>	<u>4/19</u>	Oxycodone	<u>EK 16167</u>	<u>10/18</u>

Comments: verification following maintenance.
ay 6/26/18

Batch: GIA-2060626U-AAJ

Inventory # Exp.

AMP - 17 2/19

BARB - 16 2/19

BZ - 7 3/19

BE - 4 3/19

THC - 6 3/19

CARISO - 17 4/19

PBS: 840918-PBS7 Exp. 10/9/18

Pipettes: 0395, 235D

Controls:

Blank: 101017U Exp. 10/10/18

MTDN - 16 12/18

METH - 4 11/18

OPI - 4 2/19

PCP - 20 2/19

ZOL - n/a

OXY - 7U 10/18

Mixed

Exp.

Oxy

Exp.

Neg: 050718U-MXD 10/27/18

Neg: 050718U-Q-50 10/9/18

Cal: 050418U-MXC

Cal: 050418U-C-100

Pos: 050718U-MXP

Pos: 050418U-Q-200

TMB: ~~0506126/18~~

Zol TMB: ~~0506126/18~~

STOP: _____

Notes:

Verification following maintenance. 0506126/18

0506126/18

Tuesday, June 26, 2018 10:22:01 AM

Houston

Forensic Science Center Page 1 of 2

EIA_20180626U_AAJ PBS: 040918-PBS7 Pipettes: 235D, 0395 Verification

ID	209	210	214	206	205	231	232	211	207
1 Blank-1	X	X	X	X	X	X	X	X	X
2 Blank-2	X	X	X	X	X	X	X	X	X
3 Negative-1	X	X	X	X	X	X	X	X	X
4 Negative-2	X	X	X	X	X	X	X	X	X
5 Calibrator-1	X	X	X	X	X	X	X	X	X
6 Calibrator-2	X	X	X	X	X	X	X	X	X
7 Positive-1	X	X	X	X	X	X	X	X	X
8 Positive-2	X	X	X	X	X	X	X	X	X
9 Oxy Negative-1	0	0	0	0	0	0	0	0	0
10 Oxy Negative-2	0	0	0	0	0	0	0	0	0
11 Oxy Calibrator-1	0	0	0	0	0	0	0	0	0
12 Oxy Calibrator-2	0	0	0	0	0	0	0	0	0
13 Oxy Postive-1	0	0	0	0	0	0	0	0	0
14 Oxy Postive-2	0	0	0	0	0	0	0	0	0
15 050718U-MXN	X	X	X	X	X	X	X	X	X
16 050418U-MXC	X	X	X	X	X	X	X	X	X
17 050718U-MXP	X	X	X	X	X	X	X	X	X
18 050718U-Q-50	0	0	0	0	0	0	0	0	0
19 050418U-C-100	0	0	0	0	0	0	0	0	0
20 050718U-Q-200	0	0	0	0	0	0	0	0	0

SEQUENCE VERIFIED

6/26/18

AM

AM

Tuesday, June 26, 2018 10:22:01 AM

Houston

Forensic Science Center Page 2 of 2

EIA_20180626U_AAJ PBS: 040918-PBS7 Pipettes: 235D, 0395 Verification

ID	208	221
1 Blank-1	X	X
2 Blank-2	X	X
3 Negative-1	X	O
4 Negative-2	X	O
5 Calibrator-1	X	O
6 Calibrator-2	X	O
7 Positive-1	X	O
8 Positive-2	X	O
9 Oxy Negative-1	O	X
10 Oxy Negative-2	O	X
11 Oxy Calibrator-1	O	X
12 Oxy Calibrator-2	O	X
13 Oxy Postive-1	O	X
14 Oxy Positive-2	O	X
15 050718U-MXN	X	O
16 050418U-MXC	X	O
17 050718U-MXP	X	O
18 050718U-Q-50	O	X
19 050418U-C-100	O	X
20 050718U-Q-200	O	X



QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

Exp. Group Num 2

Validation criteria : Difference data

NC2>LPC2 --> TRUE

LPC2>PC2 --> TRUE

PC2>HPC2 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8	
<i>AMP</i>	A1	209	Blank-1	3.51	3.4805	1.1987	100.85	100	
	B1	209	Blank-2	3.451			99.152		
	C1	209	Negative-1	1.813	1.728	6.9565	52.09	49.648	
	D1	209	Negative-2	1.643			47.206		
	E1	209	Calibrator-1	1.077	1.039	5.1723	30.944	29.852	
	F1	209	Calibrator-2	1.001			28.76		
	G1	209	Positive-1	0.57	0.5575	3.1709	16.377	16.018	
	H1	209	Positive-2	0.545			15.659		
	A2	209	050718U-MXN	1.528	1.528		43.902	43.902	neg
	B2	209	050418U-MXC	0.919	0.919		26.404	26.404	POS
<i>Barb</i>	C2	209	050718U-MXP	0.515	0.515		14.797	14.797	POS
	A3	210	Blank-1	3.384	3.4175	1.3863	99.02	100	
	B3	210	Blank-2	3.451			100.98		
	C3	210	Negative-1	1.626	1.5905	3.1565	47.579	46.54	
	D3	210	Negative-2	1.555			45.501		
	E3	210	Calibrator-1	1.266	1.2425	2.6748	37.045	36.357	
	F3	210	Calibrator-2	1.219			35.669		

EJA-20180626U-AAJ

clay

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
G3	210	Positive-1	0.89	0.9165	4.0891	26.042	26.818	
H3	210	Positive-2	0.943			27.593		
A4	210	050718U-MXN	1.583	1.583		46.32	46.32	neg
B4	210	050418U-MXC	1.364	1.364		39.912	39.912	neg
C4	210	050718U-MXP	0.913	0.913		26.715	26.715	POS ✓

EIA-20180626U-AAS

aaj

QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

Exp. Group Num 2

Validation criteria : Difference data

NC2>LPC2 --> TRUE

LPC2>PC2 --> TRUE

PC2>HPC2 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
BZ A1	214	Blank-1	2.797	2.749	2.4693	101.75	100	
B1	214	Blank-2	2.701			98.254		
C1	214	Negative-1	1.035	1.088	6.8891	37.65	39.578	
D1	214	Negative-2	1.141			41.506		
E1	214	Calibrator-1	0.55	0.5715	5.3203	20.007	20.789	
F1	214	Calibrator-2	0.593			21.571		
G1	214	Positive-1	0.311	0.3095	0.6854	11.313	11.259	
H1	214	Positive-2	0.308			11.204		
A2	214	050718U-MXN	1.095	1.095		39.833	39.833	neg
B2	214	050418U-MXC	0.421	0.421		15.315	15.315	POS
C2	214	050718U-MXP	0.311	0.311		11.313	11.313	POS
BE A3	206	Blank-1	2.493	2.5265	1.8752	98.674	100	
B3	206	Blank-2	2.56			101.33		
C3	206	Negative-1	1.83	1.8155	1.1295	72.432	71.858	
D3	206	Negative-2	1.801			71.284		
E3	206	Calibrator-1	1.599	1.617	1.5743	63.289	64.002	
F3	206	Calibrator-2	1.635			64.714		

EIA-20180626U-AAA5

Ray

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
G3	206	Positive-1	1.447	1.407	4.0205	57.273	55.69	
H3	206	Positive-2	1.367			54.106		
A4	206	050718U-MXN	1.85	1.85		73.224	73.224	neg
B4	206	050418U-MXC	1.565	1.565		61.943	61.943	POS
C4	206	050718U-MXP	1.442	1.442		57.075	57.075	POS

ETA-20180626U-AAJ

QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

Exp. Group Num 2

Validation criteria : Difference data

NC2>LPC2 --> TRUE

LPC2>PC2 --> TRUE

PC2>HPC2 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
<i>THE</i> A1	205	Blank-1	2.964	2.9765	0.59391	99.58	100	
B1	205	Blank-2	2.989			100.42		
C1	205	Negative-1	1.977	1.9735	0.25081	66.42	66.303	
D1	205	Negative-2	1.97			66.185		
E1	205	Calibrator-1	1.57	1.512	5.4249	52.747	50.798	
F1	205	Calibrator-2	1.454			48.849		
G1	205	Positive-1	0.802	0.8065	0.78908	26.944	27.096	
H1	205	Positive-2	0.811			27.247		
A2	205	050718U-MXN	2.109	2.109		70.855	70.855	neg
B2	205	050418U-MXC	1.605	1.605		53.922	53.922	neg
C2	205	050718U-MXP	0.786	0.786		26.407	26.407	POS
<i>Cariso</i> A3	231	Blank-1	3.096	3.1255	1.3348	99.056	100	
B3	231	Blank-2	3.155			100.94		
C3	231	Negative-1	1.98	2.0365	3.9235	63.35	65.158	
D3	231	Negative-2	2.093			66.965		
E3	231	Calibrator-1	1.662	1.658	0.34119	53.175	53.048	
F3	231	Calibrator-2	1.654			52.92		

QIA-20180626U-PAJ

QAJ

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
G3	231	Positive-1	1.214	1.201	1.5308	38.842	38.426	
H3	231	Positive-2	1.188			38.01		
A4	231	050718U-MXN	2.067	2.067		66.133	66.133	neg
B4	231	050418U-MXC	1.674	1.674		53.559	53.559	neg
C4	231	050718U-MXP	1.299	1.299		41.561	41.561	POS



GIA-20180626U-AAJ

QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

Exp. Group Num 2

Validation criteria : Difference data

NC2>LPC2 --> TRUE

LPC2>PC2 --> TRUE

PC2>HPC2 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
<i>MTDH</i> A1	232	Blank-1	1.74	✓ 1.758	✓ 1.448	98.976	100	
B1	232	Blank-2	1.776			101.02		
C1	232	Negative-1	0.495	0.492	0.86233	28.157	27.986	
D1	232	Negative-2	0.489			27.816		
E1	232	Calibrator-1	0.246	0.2565	5.7892	13.993	14.59	
F1	232	Calibrator-2	0.267			15.188		
G1	232	Positive-1	0.138	0.1395	1.5207	7.8498	7.9352	
H1	232	Positive-2	0.141			8.0205		
A2	232	050718U-MXN	0.593	0.593		33.732	33.732	neg
B2	232	050418U-MXC	0.313	0.313		17.804	17.804	neg ✓
C2	232	050718U-MXP	0.14	0.14		7.9636	7.9636	POS ✓
<i>Meth</i> A3	211	Blank-1	2.398	✓ 2.4485	✓ 2.9168	97.938	100	✓
B3	211	Blank-2	2.499			102.06		
C3	211	Negative-1	1.153	1.128	3.1343	47.09	46.069	
D3	211	Negative-2	1.103			45.048		
E3	211	Calibrator-1	0.742	0.7665	4.5203	30.304	31.305	
F3	211	Calibrator-2	0.791			32.305		

EIA-20180626U-AAJ

AAJ

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
G3	211	Positive-1	0.565	0.5675	0.623	23.075	23.177	
H3	211	Positive-2	0.57			23.28		
A4	211	050718U-MXN	0.918	0.918		37.492	37.492	neg
B4	211	050418U-MXC	0.721	0.721		29.447	29.447	POS
C4	211	050718U-MXP	0.434	0.434		17.725	17.725	POS



ETA-20180626U-PAJ

QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

Exp. Group Num 2

Validation criteria : Difference data

NC2>LPC2 --> TRUE

LPC2>PC2 --> TRUE

PC2>HPC2 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
<i>OPI</i> A1	207	Blank-1	3.099	3.0925	0.29725	100.21	100	
B1	207	Blank-2	3.086			99.79		
C1	207	Negative-1	0.514	0.5485	8.8952	16.621	17.736	
D1	207	Negative-2	0.583			18.852		
E1	207	Calibrator-1	0.317	0.329	5.1582	10.251	10.639	
F1	207	Calibrator-2	0.341			11.027		
G1	207	Positive-1	0.199	0.189	7.4826	6.4349	6.1116	
H1	207	Positive-2	0.179			5.7882		
A2	207	050718U-MXN	0.602	0.602		19.466	19.466	neg
B2	207	050418U-MXC	0.309	0.309		9.9919	9.9919	POS
C2	207	050718U-MXP	0.176	0.176		5.6912	5.6912	POS
<i>PCP</i> A3	208	Blank-1	2.951	3.0025	2.4257	98.285	100	
B3	208	Blank-2	3.054			101.72		
C3	208	Negative-1	1.384	1.406	2.2129	46.095	46.828	
D3	208	Negative-2	1.428			47.56		
E3	208	Calibrator-1	0.906	0.9155	1.4675	30.175	30.491	
F3	208	Calibrator-2	0.925			30.808		

EIA-20180626-011-AAJ

AAJ

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
G3	208	Positive-1	0.619	0.6095	2.2043	20.616	20.3	
H3	208	Positive-2	0.6			19.983		
A4	208	050718U-MXN	1.474	1.474		49.092	49.092	neg
B4	208	050418U-MXC	0.943	0.943		31.407	31.407	neg
C4	208	050718U-MXP	0.575	0.575		19.151	19.151	POS ✓

CTA-20180626U-AAJ

RAJ

QC Validation criteria

Exp. Group Num 1

QC Validation criteria : Difference data

NC1>LPC1 --> TRUE

LPC1>PC1 --> TRUE

PC1>HPC1 --> TRUE

1 Strip method names

2 Sample ID 1

3 Difference data

4 Difference data - Mean

5 Difference data - Variation coefficient

6 b/b0

7 b/b0 - Mean

8 Cutoff results

	1	2	3	4	5	6	7	8
OXY A1	221	Blank-1	3.278	✓ 3.2635	✓ 0.62835	100.44	100	
B1	221	Blank-2	3.249			99.556		
C1	221	Oxy Negative-1	1.519	1.513	0.56082	46.545	46.361	
D1	221	Oxy Negative-2	1.507			46.177		
E1	221	Oxy Calibrator-1	1.34	1.3355	0.47652	41.06	40.922	
F1	221	Oxy Calibrator-2	1.331			40.784		
G1	221	Oxy Postive-1	1.064	1.0605	0.46674	32.603	32.496	
H1	221	Oxy Positive-2	1.057			32.389		
A2	221	050718U-Q-50	1.564	1.564		47.924	47.924	neg
B2	221	050418U-C-100	1.332	1.332		40.815	40.815	POS
C2	221	050718U-Q-200	1.08	1.08		33.093	33.093	POS

EIA-20180626U-AAJ

AAJ