

**HOUSTON POLICE DEPARTMENT**  
**CRIME LABORATORY**  
*Toxicology Section*

**Reagent Preparation Worksheet**  
**Volatiles Analysis**

**Reagent:** 1.0% n-propanol Stock Solution

**Preparation Notes:** In a 100 ml volumetric flask, weight out 1 g of n-propyl alcohol (n-propanol) and bring to volume with DI water

Lab Lot #	Date Prepared	Preparer's Initials	Chemical 1		Chemical 2		Chemical 3		Volume		QC Check		Expiration Date	
			Manufacturer/Lot #	Weight	Manufacturer/Lot #	Weight	Manufacturer/Lot #	Weight	Made	Date	Date	Date		
01102014	1/16/14	DW	n-propanol Sigma-Aldrich Lot#: 503910 LMV	weighed out 1.0024 g	Millipore Direct Q W3 DI H2O				100ml		1/21/14		7/10/14	