



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

Working Stock/Standard Preparation Log

Solution Name/Concentration: Oxycodone Blood Cut-off Calibrator / 10ng/mL

Solvent/Matrix Used: Blood Immalysis E34447
Solvent/Matrix Name Manufacturer Lot Number

New Solution: 180924-B-C-10 ay 9/24/2018
Lot Number Prepared By Date

3/19/2019 CIA-20180926B-AAJ
Expiration Date Verification Batch Name ay 9/28/18

Pipette(s): 5581

Balance(s): Analytical: 1123241782 Top Loading: 1123253484 (N/A)

Volumetric Flask SN#: N/A

Preparation Notes: Add 50 uL of Oxycodone cut-off calibrator stock to a 25 mL volumetric flask. DS with blood and thoroughly mix.

#	Drug Standard/Concentration	Expiration Date	Inventory Number	Manufacturer	Amount Added	Lot Number	Final Concentration
1	Oxycodone cut-off Calibrator Stock / 5 ng/mL	3/19/19	N/A	In-house	50 uL	180924C-5-5	10 ng/mL
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							