



# Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

## Working Stock/Standard Preparation Log

Solution Name/Concentration: Oxycodone Cut-off Calibrator Stock / 5ug/mL

Solvent/Matrix Used: PBS In-house 180919-PBS7  
Solvent/Matrix Name Manufacturer Lot Number

New Solution: 180924C-SE-5 ag 9/24/18  
Lot Number Prepared By Date

3/19/2019 In dilution 180924B-C-10  
Expiration Date Verification Batch Name orig 9/24/18

Pipette(s): 5581

Balance(s): Analytical: 1123241782 Top Loading: 1123253484 (N/A)

Volumetric Flask SN#: N/A

Preparation Notes: Add 50 uL of Oxycodone to a 10mL volumetric flask and QS with PBS.

#	Drug Standard/Concentration	Expiration Date	Inventory Number	Manufacturer	Amount Added	Lot Number	Final Concentration
1	Oxycodone / 1ug/mL	6/2021	17-43	Cerilliant	50uL	FE06081603	5ug/mL
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							

orig 9/24/18