



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: GCMS-2

Date	He Tank Pressure	H ₂ Tank Pressure*	Air Tank Pressure*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Bead Voltage*	Comments	Signature
11/19/18	70 2200	X	X	X	X	X	X	X	✓	✓	✓	n/a	n/a	clipped about a foot of the column. NO casework was run today.	TS
11/20/18	70 2200	X	X	✓	X	✓	✓	✓	X	X	X	n/a	n/a		TS
mk 11/20/18															
 															
 															
 															

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melina Henney

Date Completed: 11/20/18

