



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: 11/26/18 mh
#2 GCMS 2

Date	He Tank Pressure	H ₂ Tank Pressure*	Air Tank Pressure*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Bead Voltage*	Comments	Signature
11/20/18	n/a	n/a	n/a	✓	X	X	✓	✓	✓	✓	X	n/a	n/a	Added PFTBA to cal vial. Changed rough pump oil + diff pump fluid. Cleaned MS fans. Cleaned source. Changed inlet liner + O-ring. Changed gold seal and washer. Cleaned split vent copper tubing. Cleaned MS fan. use 11-26-18 cleaned ALS guide rods + rotated	
														betts. Cleaned tray guide rods.	L. Busch FSA mh 11/26/18
11/27/18	70 1700	n/a	n/a	✓	X	X	✓	✓	X	X	X	n/a	n/a	post PM autotune	mh

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melvin Henry

Date Completed: 12/3/18

