



**Houston Forensic Science Center**  
Comparative and Analytical Division - Toxicology

**GC/MS Maintenance Log**

Instrument: GCMS-3

Date	He Tank Pressure	Air Tank Pressure*	H <sub>2</sub> Water Level*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Comments	Analyst
			H <sub>2</sub> Pressure*									Bead Voltage*		
08/31/2020	70 1500	nta	nta	X	X	X	X	✓	X	X	X	nta	FSA rep onsite; MS for line pressure normal and dropping. Following up tomorrow.	DMC
09/01/2020	70 1500	nta	nta	✓	X	X	X	✓	X	X	X	nta		DMC
9/3/2020 m/h														

\*For GC/MS-NPD instrument (GCMS-5).

Signature: Melanie Henry

Date Completed: 9/3/2020

