



Houston Forensic Science Center
Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: GCMS-4

Date	He Tank Pressure	H ₂ Tank Pressure*	Air Tank Pressure*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Bead Voltage*	Comments	Signature
9/18/18	/	/	/	✓										Lisa Busch from Full spectrum changed the copper tubing, split vent filter, furren on the column.	TS
9/19/18	70 2300	n/a	n/a	✓	X	X	X	✓	X	X	X	n/a	n/a		TS
9/20/18	70 2250	n/a	n/a	✓	X	✓	✓	✓	X	X	X	n/a	n/a		mh
9/21/18	70 2200	n/a	n/a	✓	✓	✓	✓	✓	✓	X	X	n/a	n/a		TS
9/28/18 mb															

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melvin Henry

Date Completed: 9/28/18

