



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: GCMS-4

Date	He Tank Pressure	H ₂ Tank Pressure*	Air Tank Pressure*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Bead Voltage*	Comments	Signature
9/24/18	70 2100	n/a	n/a	✓	✓	✓	✓	✓	✓	X	X	n/a	n/a		TS
9/26/18	70 1900	n/a	n/a	X n/a 9/26/18 mh	X	X	X	X	✓	✓	✓	n/a	n/a	clipped columns cleaned source	mh
9/27/18	70 1800	n/a	n/a	✓	✓	✓	✓	✓	X	X	X	n/a	n/a	replaced syringe	mh
9/28/18	70 1700	n/a	n/a	✓	X	✓	✓	✓	X	X	X		n/a		mh
9/28/18 mh															

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melina Henry

Date Completed: 9/28/18

LAB-24

Issued By: Supervisor - Toxicology

Issue Date: 2017-03-06

