



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: GCMS-4

Date	He Tank Pressure	H ₂ Tank Pressure*	Air Tank Pressure*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Bead Voltage*	Comments	Signature
11/13/18	70 700	X	X	✓	✓	✓	✓	✓	✓	X	X	N/A	N/A		luc
11/14/18	70 600	X	X	✓	X	✓	✓	✓	✓	X	X	N/A	N/A		luc
11/15/18	70 2500	X	X	✓	✓	✓	✓	✓	✓	X	X	N/A	N/A	Chemistis	luc
11/16/18	70 2500	X	X	✓ * TS 11/16/18	✓	X	X	X	✓	✓	✓	Na	n/a	about a foot of the column was clipped & nothing was run today. No casework was run today.	TS
11/16/18 mh															

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melina Henry

Date Completed: 11/16/18

LAB-24

Issued By: Supervisor - Toxicology

Issue Date: 2017-03-06

