



Houston Forensic Science Center

Comparative and Analytical Division - Toxicology

GC/MS Maintenance Log

Instrument: GCMS-5

Date	He Tank Pressure	Air Tank Pressure*	H ₂ Water Level*	Auto Tune	Septum Change	Wash Vials Filled	Syringe Washed	Rough Pump Oil Level Checked	Liner Change	Gold Seal Replaced	Column Cut	Output*	Comments	Analyst
			H ₂ Pressure*									Bead Voltage*		
8/25/2020													Shutting down in preparation for hurricane. -VC 8/25/2020	VC
8/27/2020	90 1350	80 700	yes 80	✓	X	X	X	✓	X	X	X	0.4 0.6B	started back up. baked out MS.	mlh
8/28/2020	90 1350	80 700	yes 80	✓	X	X	X	✓	X	X	X	0.4 0.6B		mlh
8/28/2020 mlh														

*For GC/MS-NPD instrument (GCMS-5).

Signature: Melvin Henry

Date Completed: 8/28/2020

